
# WILOVERT

INVERSE MIKROSKOP-REIHE

INVERTED MICROSCOPES

## INVERSE MIKROSKOPE

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## INVERTED MICROSCOPES

have gradually become an indispensable tool in many laboratories where living material is examined. Also, they are being increasingly upgraded into highly specialized systems for certain routine applications and particular research developments. The new Wilovert microscopes are especially designed for this trend, featuring:

### HIGH STABILITY AND ERGONOMY

The entire U-shaped body of the Wilovert consists of one cast part, which provides optimal conditions for high stability. Micromanipulators can easily be adapted to any microscope in the Wilovert series.

### WIDE CHOICE OF LIGHT SOURCES

Besides the 6V / 20W and the 6V / 30 W halogen lamps there is also a choice of 12V / 50W or 12V/100W lamps which can also be adapted from the outside. For fluorescence work we offer HBO 50W or HBO 100W mercury lamps.

### OPTICS

**Objectives**  
As inverted microscope specimens are examined in tissue culture vessels of different heights, base thicknesses and materials, as opposed to the standard specimen slides and coverglasses used for upright microscopes, special optics are required. Hund offers objectives and matching condensers with long free working distances.

### CONDENSERS

There are the following condensers to choose from:

- Condenser NA 0.25 for transmitted light brightfield with 58 mm free working distance.
- Combined condenser NA 0.25 with 58 mm free working distance for transmitted light brightfield and phase contrast.
- Condenser NA 0.4 for transmitted light brightfield with 55 mm free working distance.
- Combined condenser NA 0.4 with 38 mm free working distance for transmitted light brightfield and phase contrast.

### HOLDERS FOR CULTURE VESSELS

Frequently used culture vessels such as petri dishes, tissue culture vessels, microtest trays and microtiter trays are inserted in special holding frames on Wilovert microscopes. The object guide ensures that specimens are completely scanned.

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### WILOVERT AS HF

The standard outfit for routine examinations of transparent specimens with richly contrasting details which can be done in brightfield. The advantage of this outfit is the simple, clearly arranged configuration which leaves the user plenty of room for manoeuvre.

#### CONFIGURATION WILOVERT AS HF:

Brightfield objectives	A 4/0.10
	A 10/0.25
Pair of eyepieces	WF 10x/18

### WILOVERT AS PH

Routine microscope for monitoring cell cultures in brightfield and phase contrast.

#### CONFIGURATION WILOVERT AS PH:

Brightfield objective	A 4/0.10
Brightfield/PH objective	A Ph 10/0.25
Brightfield/PH objective	SPL Ph 20/0.35
Pair of eyepieces	WF 10x/18
Combined condenser NA 0.25 (brightfield + phase contrast)	

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### WILOVERT SE HF

The standard configuration for examinations in transmitted light brightfield

#### CONFIGURATION WILOVERT SE HF:

Brightfield objectives	A 4/0.10
	A 10/0.25
	SPL 20/0.35

Condenser NA 0.25, pair of eyepieces WF 10x/18  
Object guide with left-hand drive

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### WILOVERT SE PH

The standard microscope for switching between brightfield or phase contrast.

#### CONFIGURATION WILOVERT SE PH:

Same optics as Wilovert SE HF, but with objectives A Ph 10x, SPL Ph 20/0.35 with combined condenser NA 0.25 (brightfield and phase contrast) and focusing telescope. (see illustration on front cover)

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### WILOVERT S LL 50/100

Illumination is a particularly important component of inverted microscopes, as it has to satisfy such a wide variety of demands.

#### **HERE, HUND MAKES USE OF ITS COMPETENCE IN FIBRE OPTICS TECHNOLOGY.**

Using a 12 V/50 W or 12 V/100 W halogen lamp in an external lamphousing, light is directed to the illumination system via fibre optics. This has the following advantages:

- The specimen is not exposed to any heat.
- There are no electric or magnetic fields to influence e.g. electrophysiological examinations.
- There is no heavy lamphousing weighing down the transmitted light illumination column.

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### WILOVERT S HF-LL

Outfit for scientific work and highly qualified routine examinations in brightfield.

#### **CONFIGURATION WILOVERT S HF-LL:**

Brightfield objectives	A	4/0.10
	A	10/0.25
	SPL	20/0.35

Condenser NA 0.25, pair of eyepieces WF 10x/18

Object guide with left-hand drive

With external lamphousing HAL 50, optionally available with lamphousing HAL 100.

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### WILOVERT S PH-LL

Outfit for scientific work and highly qualified routine examinations in brightfield and phase contrast.

#### **CONFIGURATION WILOVERT S PH-LL:**

Brightfield objective	A	4/0.10
Brightfield/PH objective	A PH	10/0.25
Brightfield/PH objective	SPL PH	20/0.35

Combined condenser NA 0.25 (brightfield + phase contrast), focusing telescope, pair of eyepieces WF 10x/18, left-hand object guide.

With external lamphousing HAL 50, optionally available with lamphousing HAL 100.

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## WILOVERT AFL FOR EXAMINATIONS IN INCIDENT LIGHT FLUORESCENCE


Fluorescence is a key technology in the development of molecular-biological and biochemical examination and measurement techniques, e.g. for visualising biomolecules in their natural environment (in cells and tissue) or for studying physiological processes. To obtain bright, high-contrast fluorescence images, the Wilovert AFL features the ultra high pressure mercury lamps (HBO 50, HBO 100) generally preferred for fluorescence microscopy because of their high performance. Filter cubes are individually interchangeable, and any of those in the Hund product range can be chosen.

## WILOVERT S HF AFL FOR OBSERVATIONS IN TRANSMITTED LIGHT AND INCIDENT LIGHT FLUORESCENCE, E.G. FOR HLA TYPING

### CONFIGURATION WILOVERT S HF AFL:

Brightfield objectives A 4/0.10  
A 10/0.25  
SPL 20/0.35

Condenser NA 0.25  
Pair of eyepieces WF 10x/18  
Object guide with left-hand drive  
With external lamphousing HBO 50, optionally available with lamphousing HBO 100.

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## WILOVERT S PH AFL UNIVERSALLY APPLICABLE OUTFIT FOR OBSERVATIONS IN INCIDENT LIGHT FLUORESCENCE IN COMBINATION WITH TRANSMITTED LIGHT PHASE CONTRAST, E.G. FOR ENABLING CLASSIFICATION OF FLUORESCING STRUCTURES IN THEIR ENVIRONMENT.

### CONFIGURATION WILOVERT S PH AFL:

Brightfield objective A 4/0.10  
Brightfield/PH objective A PH 10/0.25  
Brightfield/PH objective SPL PH 20/0.35  
Combined condenser NA 0.25 (brightfield and phase contrast), focusing telescope  
Pair of eyepieces WF 10x/18  
Object guide with left-hand drive  
With external lamphousing HBO 50, optionally available with lamphousing HBO 100.

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## ACCESSORIES

### Micromanipulators and microinjectors

Micromanipulators and microinjection systems are now commercially available for working on living cells – e.g. for injecting substances into cell nuclei or cytoplasm. The stage of Wilovert microscopes proves ideal for work with these accessories. Due to the importance of being vibration-free for microinjection and electrophysiological tests, it is fixed on both sides of the microscope stand.

As focusing is performed by sliding the objective nosepiece, not by raising or lowering the stage, the system stage-specimen-micromanipulator (microinjector) always maintains a constantly high level of stability.

## UPGRADE OPTIONS IMAGE DOCUMENTATION

The following accessories can be used for documenting microscope images:

- Photographic systems
- Video systems
- Image archiving

There is a choice of two different binocular phototubes for adapting the documentary equipment via standard components.

## BINOCULAR PHOTOTUBES

Binocular phototube 70/30

70% of the light is directed to the photographic or video camera, 30% to the eyepieces for simultaneous viewing.

Binocular phototube 100/100

The light can either be fully directed to the eyepieces or to the connected documentary equipment. This is a particular advantage for dark images as the full light intensity is available. Video cameras are connected via C-mount adapter.

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Helmut Hund GmbH